
mctOptics Documentation

Release 0.1

Argonne National Laboratory

Apr 04, 2023

CONTENTS


1	Content	1
----------	----------------	----------

CHAPTER ONE

CONTENT

1.1 About

mctOptics is an EPICS IOC supporting the Optique Peter system installed at beamline 2-BM of the Advanced Photon Source.




1.2 Instrument

1.2.1 Optique Peter

Triple Objective

At beamline 2-BM we use a Triple Objective Microscope produced by [Optique Peter](#)



Detailed information of the instruments are:

1. [Test Report](#)
2. [User manual](#)
3. [Manuals](#)
4. [Specs](#)
5. [Reference documentation](#)

Lens	Lens Mag	motor position (specs)	ref on lens 1 (specs)	ref on lens 1 (aligned)
0	10 x	121.5942	59.6099	59.0151
1	5 x	61.9841	0.0000	-0.3690
2	1.1 x	2.3006	-59.6835	-59.6865

With FLIR Oryx ORX-10G-51S5M-C


Lens	Lens Mag	Pixel size (m) (specs)	Pixel size (m) (meas.)	Pixels (H)	Pixels (V)	Detector pixel size (m)	H filed of view (mm)	V filed of view (mm)
0	10 x	0.35	0.350787	2448	2048	3.45	0.84	0.71
1	5 x	0.69	0.699447	2448	2048	3.45	1.69	1.41
2	1.1 x	3.14	3.125830	2448	2048	3.45	7.68	6.42

Scintillators:

Lens	Lens Mag	Scintillator material	Scintillator thickness (m)
0	10 x	LuAG	25
1	5 x	LuAG	50
2	1.1 x	LuAG	100

1.3 Usage

Objective and camera change can be accomplished by simply selecting the desired magnification and the camera in the user interface selector of the main **mctOptics** control screen:



When changing lens/camera, **mctOptics** is also correcting for minor miss-alignment of the instrument visible light optics so that sample point of interest stays in the center of the image at each lens change. **mctOptics** also keeps the

rotation axis aligned with the detector columns and each lens/camera change by rotating the camera.

This capability allows for step-zoom-in during a tomographic measurement and shown in [this video](#)

The required lens offset sample x, y, z and the lens offset camera rotation are very reproducible and can be determined once when the instrument is first installed.

1.3.1 Energy change

The beamline x-ray energy change is managed by the `energy cli` python library.


For user operation the `energy cli` is called from the `mctOptics` user interface:



The DMM energy change operates in two modes. The first uses pre-stored energy calibration files. To select this mode set to “No” the “Use File” button:



Then you can select any available energy from the drop down list:



Once the desired energy is selected press the “Set” button to move the DMM.

The second mode allows the use of arbitrary configuration files. You can enable this mode by selecting “Yes” in the “Use File” button:



In this example we are setting the DMM to 17.04 keV using the `energy2bm_Mono_17.04.conf` configuration file. Press the “Set” button to move the DMM to this energy.

You can create new configuration files with the `energy` cli by moving the DMM to a known energy, i.e. using calibration foils, and storing that DMM positions in a new config file with:

```
[user2bmb@arcturus,42,~]$ energy save --mode Mono --energy-value 27.00
```

in this case a new config file called `energy2bm_Mono_27.0_2022-11-03_23_26_17.conf` is automatically generated and it can be used by updating the “Folder” and “File” entries in the user interface.

Finally, if you want to set the DMM to an arbitrary energy that is not a calibration point, let’s say 27.18, you can use the “interpolate” function. To do so, enter an arbitrary energy within the energy calibrated range, in our case between 11.10 and 33.20 keV, in the user interface.

The interpolation function will create a config file for 27.18 keV by interpolating linearly all DMM positions between the 2 closest calibrated energies. For 27.18 keV, the interpolation will occur using 25.51 and 30.00 keV. Once the new `energy2bm_interp_27.18.conf` is generated, its name will be automatically copied in the “File” entry of the user interface. To move the DMM to this energy press the “Set” button.



More information on how to operate the mctOptics user interface is [here](#).

1.4 Demo

mctOptics allows for step-zoom-in during a tomographic measurement and shown in [this video](#).

1.5 Install directions

1.5.1 Build EPICS base

Warning: Make sure the disk partition hosting ~/epics is not larger than 2 TB. See [tech talk](#) and [Diamond Data Storage](#) document.

```
$ mkdir ~/epics  
$ cd epics
```

- Download EPICS base latest release, i.e. 7.0.3.1., from <https://github.com/epics-base/epics-base>:

```
$ git clone https://github.com/epics-base/epics-base.git  
$ cd epics-base  
$ make -sj
```

1.5.2 Build a minimal synApps

To build a minimal synApp:

```
$ cd ~/epics
```

- Download in ~/epics `assemble_synApps.sh`
- **Edit the `assemble_synApps.sh` script as follows:**
 1. Set `FULL_CLONE=True`
 2. Set `EPICS_BASE` to point to the location of EPICS base. This could be on APSshare (the default), or a local version you built.

For mctoptics you need

1. ASYN=R4-37
2. AUTOSAVE=R5-10
3. BUSY=R1-7-2
4. XXX=R6-1

You can comment out all of the other modules (ALLENBRADLEY, ALIVE, etc.)

- Run:

```
$ assemble_synApps.sh
```

- This will create a synApps/support directory:

```
$ cd synApps/support/
```

- Edit asyn-RX-YY/configure/RELEASE to comment out the lines starting with:

```
IPAC=$(SUPPORT)/
SNCSEQ=$(SUPPORT)/
```

Warning: If building for RedHat8 uncomment **TIRPC=YES** in asyn-RX-YY/configure/CONFIG_SITE

- Clone the mctoptics module into synApps/support:

```
$ git clone https://github.com/tomography/mctoptics.git
```

- Edit configure/RELEASE add this line to the end:

```
MCTOPTICS=$(SUPPORT)/mctoptics
```

- Edit Makefile add this line to the end of the MODULE_LIST:

```
MODULE_LIST += MCTOPTICS
```

- Run the following commands:

```
$ make release
$ make -sj
```

1.5.3 Testing the installation

- Edit /epics/synApps/support/mctoptics/configure to set EPICS_BASE to point to the location of EPICS base, i.e.:

```
EPICS_BASE=/APSshare/epics/base-3.15.6
```

- Start the epics ioc and associated medm screen with:

```
$ cd ~/epics/synApps/support/mctoptics/iocBoot/iocMCTOptics
$ start_IOC
$ start_medm
```

1.6 mctOpticsApp EPICS application

mctOptics includes a complete example EPICS application, including:

- A database file and corresponding autosave request file that contain the PVs required by the mctoptics.py base class.
- OPI screens for medm
- An example IOC application that can be used to run the above databases. The databases are loaded in the IOC with the example substitutions file, mctOptics.substitutions.

1.6.1 Base class files

The following tables list all of the records in the mctOptics.template file. These records are used by the mctoptics base class and so are required.

mctOptics.template

This is the database file that contains only the PVs required by the mctoptics.py base class mctOptics.template.

TomoScan and Camera PV Prefixes

Record Name	Description
\$(P)\$(R)Detector0PVPrefix	prefix for the detector 0, e.g. 2bmbSP1:
\$(P)\$(R)Detector1PVPrefix	prefix for the detector 1, e.g. 2bmbSP2:
\$(P)\$(R)OverlayPlugin0PVPrefix	prefix for OverlayPlugin 0, e.g. 2bmbSP1:Over1:
\$(P)\$(R)OverlayPlugin1PVPrefix	prefix for OverlayPlugin 1, e.g. 2bmbSP2:Over1:
\$(P)\$(R)FilePlugin0PVPrefix	prefix for FilePlugin 0, e.g. 2bmbSP1:HDF1:
\$(P)\$(R)FilePlugin1PVPrefix	prefix for FilePlugin 1, e.g. 2bmbSP2:HDF1:

Lens Sample X-Y-Z PV Names

Record Name	Description
\$(P)\$(R)LensSampleXPVName	for LensSampleX , e.g. 2bmS1:m2
\$(P)\$(R)LensSampleYPVName	for LensSampleY, e.g. 2bmb:25
\$(P)\$(R)LensSampleZPVName	for LensSampleZ, e.g. 2bmS1:m1

Lens Focus PV Names

Record	Description
name	type
\$(P)\${(R)lensFocus}	PVName prefix for Lens0Focus, e.g. 2bmb:m2
\$(P)\${(R)lensFocus}	PVName prefix for Lens1FocusPVName, e.g. 2bmb:m3
\$(P)\${(R)lensFocus}	PVName prefix for Lens2FocusPVName, e.g. 2bmb:m4

Camera rotation PV Names

Record	Description
name	type
\$(P)\${(R)cameraRotation}	PVName for Camera0Rotation , e.g. 2bmb:m7
\$(P)\${(R)cameraRotation}	PVName for Camera1Rotation , e.g. 2bmb:m8

Optique Peter camera selector

Record	Description
name	type
\$(P)\${(R)cameraSelector}	selector for Pos0 and Pos1 position
\$(P)\${(R)cameraSelector}	status for Camera0 and Camera1 position
\$(P)\${(R)cameraPosition}	for the Camera0
\$(P)\${(R)cameraPosition}	for the Camera1
\$(P)\${(R)cameraName}	for Pos0, e.g. Adimec
\$(P)\${(R)cameraName}	for Pos1, e.g. Flir
\$(P)\${(R)cameraMotorPVName}	motor PV name, e.g. 2bmb:m5

Optique Peter camera rotation

Record	Description
name	type
\$(P)\${(R)camera0Rotation}	Lens 0 rotation value
\$(P)\${(R)camera1Rotation}	Lens 1 rotation value
\$(P)\${(R)camera1Rotation}	Lens 1 rotation value
\$(P)\${(R)camera0Rotation}	Lens 0 rotation value
\$(P)\${(R)camera1Rotation}	Lens 1 rotation value
\$(P)\${(R)camera2Rotation}	Lens 2 rotation value

Optique Peter lens focus

Record	Description
nametype	
\$(P)\$(R)CaMVeratdLns0	Camera 0 Lens 0 focus value
\$(P)\$(R)CaMVeratdLns0	Camera 0 Lens 1 focus value
\$(P)\$(R)CaMVeratdLns0	Camera 0 Lens 2 focus value
\$(P)\$(R)CaMVeratdLns1	Camera 1 Lens 0 focus value
\$(P)\$(R)CaMVeratdLns1	Camera 1 Lens 1 focus value
\$(P)\$(R)CaMVeratdLns1	Camera 1 Lens 2 focus value

Optique Peter lens selector

Record	Description
nametype	
\$(P)\$(R)CaMVerSel	selector for Pos0 and Pos1 position
\$(P)\$(R)CaMVerSel	for the first lens
\$(P)\$(R)CaMVerSel	for the second lens
\$(P)\$(R)CaMVerSel	for the third lens
\$(P)\$(R)LensNarhSel	for Pos0, e.g. Lens0
\$(P)\$(R)LensNarhSel	for Pos1, e.g. Lens1
\$(P)\$(R)LensNarhSel	for Pos2, e.g. lens2
\$(P)\$(R)CaMVerMotPMName	motor PV name, e.g. 2bmb:m1

Detector image cross

Record	Description
nametype	
\$(P)\$(R)CaCrossSelect	

Optique Peter lens 1 offsets

Record	Description
nametype	
\$(P)\$(R)Camera0Lens1XOffset	
\$(P)\$(R)Camera0Lens1YOffset	
\$(P)\$(R)Camera0Lens1ZOffset	
\$(P)\$(R)Camera1Lens1XOffset	
\$(P)\$(R)Camera1Lens1YOffset	
\$(P)\$(R)Camera1Lens1ZOffset	

Optique Peter lens 2 offsets

Record	Record	Description
name	type	
\$(P)	\$(R)	Camera0Lens2XOffset
\$(P)	\$(R)	Camera0Lens2YOffset
\$(P)	\$(R)	Camera0Lens2ZOffset
\$(P)	\$(R)	Camera1Lens2XOffset
\$(P)	\$(R)	Camera1Lens2YOffset
\$(P)	\$(R)	Camera1Lens2ZOffset

MCT status via Channel Access

Record	Record	Description
name	type	
\$(P)	\$(R)	MCTStatus
	form	
\$(P)	\$(R)	Watchdog
	cout	
\$(P)	\$(R)	ServerRunning

Sync to motor

Record	Record	Description
name	type	
\$(P)	\$(R)	Sync

Optics information

Record	Record	Description
name	type	
\$(P)	\$(R)	ScintillatorType
		type of scintillator being used.
\$(P)	\$(R)	ScintillatorThickness
		thickness of the scintillator in microns.
\$(P)	\$(R)	ImagePixelSize
		pixel size on the sample in microns (i.e. includes objective magnification)
\$(P)	\$(R)	DetectorPixelSize
		pixel size of the detector.
\$(P)	\$(R)	CameraObjective
		objective of the camera objective
\$(P)	\$(R)	DetectorObjective
		objective of the camera objective

Lens name

Record	Description
name	type
\$(P)\${(R)lensName}	Name for Lens0, e.g. 1.1x
\$(P)\${(R)lensName}	Name for Lens1, e.g. 5x
\$(P)\${(R)lensName}	Name for Lens2, e.g. 10x

Camera names

Record	Description
name	type
\$(P)\${(R)camera0Name}	
\$(P)\${(R)camera1Name}	


Detector cropping

Record	Description
name	type
\$(P)\${(R)CutLeft}	
gout	
\$(P)\${(R)CutRight}	
gout	
\$(P)\${(R)CutTop}	
gout	
\$(P)\${(R)CutBottom}	
gout	
\$(P)\${(R)Cut}	
\$(P)\${(R)SuggestedAngles}	
\$(P)\${(R)SuggestedAngleStep}	

medm files

mctOptics.adl

The following is the MEDM screen `mctOptics.adl` during a scan. The status information is updating.



mctOpticsEPICS_PVs.adl

The following is the MEDM screen `mctOpticsEPICS_PVs.adl`.

If these PVs are changed tomoscan must be restarted.

